



Medicine pedestal drawer MH1

For safe storage of medicines/narcotics. Robur's medicine pedestal drawer MH1 is designed for the safe handling of medicines in confined spaces. Like the medicine rooms in hospital wards or smaller care units. The MH1 can of course be used in several different environments and needs. It can be connected to an unlimited number of units that are controlled with the same touch screen. One included function is reading event logs, i.e. who opened which drawer, size of withdrawal, day and time. But also functions such as stock balance and inventory are included as well as others. The login takes place via an RFID reader. The system has several authorization levels such as daily users, administrators and manager.



Elegant design in silver metallic paint.



Each drawer is equipped with a lock with high security. There is an RFID login to get access to the system and the drawers. The drawers have a mechanical key lock emergency opening.



Can be used as a stand alone solution or be part of a network solution.



The pedestal drawer can be equipped with up to 28 drawers in different sizes in many combinations. The pedestal drawer can be mounted on a frame at working height.

TECHNICAL SPECIFICATIONS

External dimensions

Height: 522mm, Width: 550mm, Depth: 410mm

User

The daily user takes out medicine by first identifying himself. The user chooses the medicament to be taken out and the system opens the correct drawer automatically. The user can also open the drawers with a hand scanner. Each withdrawal is registered by the system.

Administrators

The administrators sets which authorizations are to be registered on medicine and user level. The administrator registers the medicaments that shall be stored in the drawers. Furthermore, a number of parameters are specified for the medicaments such as order points, package units, sizes etc. There are also functions that handle the refill of medicaments. The system shows which drawer the medicament is stored in and the expiry date for the existing medicament.

Stand alone and/or network

The medicine pedestal drawer can be used as a stand-alone unit or in a larger network with one central administration. Even on the stand alone option, several units can be connected and run together. The pedestal drawers are then administered directly through the computer(touch screen) located in the space in question.

It is also possible to administer a larger number of pedestal drawers that are scattered in different premises, so-called network solution.

Security

The pedestal drawer is a solid steel construction with strong and secure locks on each drawer.

Unauthorized persons cannot access the content without inflicting great violence.

In the event of a power cut, for example, there is also one emergency opening of the drawers, which is then opened via a key lock.

Software

The software that runs the MH1 is developed by Robur. The software is a Windows-based solution and Robur works continuously with updates. Robur also handles updates from Microsoft Windows.

This is managed through a license agreement between the customer and Robur.

API

The software contains an API solution, which allows the system to be connected with an existing ordering system or various databases.

Language

Today Swedish and English

We reserve the right to alter the specification without prior notice.



Head office:
Robur Safe AB
Lundbyvägen 8
311 68, Slöinge
Sweden

☎ +46 346 260 260
✉ info@robursafe.com
🌐 www.robursafe.com